



**MAGNETIC RESONANCE IMAGING**


Edge your competition



**ContextVision**

# Image enhancement for Magnetic Resonance Imaging

---



GOPView® MRI2<sup>Plus</sup> provides great soft tissue contrast as well as superior noise reduction and edge sharpness.




MRI

**GOPView® MRI2<sup>Plus</sup>**

Setting new standards in reliability and image quality

- Intelligent noise suppression with preserved clinical details.
- Clear soft tissue differentiation (e.g. gray and white matter) and visualization of subtle structures.
- Superior contrast and edge enhancement.
- Advanced algorithms that correct for varying intensity levels in the image caused by inhomogeneities in the magnetic field.
- Enables higher diagnostic confidence and increased patient throughput.



The brain scan processed with GOPView® MRI2<sup>Plus</sup> shows excellent noise reduction and contrast enhancement.

## Customizable solutions and professional support

GOPView MRI2<sup>Plus</sup> is designed for seamless integration. The image processing library is a .dll file delivered in an SDK, together with parameter files adapted to optimize image quality and performance.

With a highly configurable design, GOPView MRI2<sup>Plus</sup> is optimized for the unique requirements of each application, sequence type and protocol.

Customizable for different preferences.

## Tuning tool

Gives the MRI manufacturer the possibility to optimize image quality by creating new parameter files.

The ContextVision team will provide professional support throughout the entire lifecycle.

# At your service. Every day. All over the world.



ContextVision is a medical technology company specializing in image analysis and artificial intelligence. Its cutting-edge technology helps doctors accurately interpret medical images, a crucial foundation for better diagnosis and treatment. As an industry pioneer for more than 30 years, ContextVision is taking a lead position within deep learning, the latest artificial intelligence technology. ContextVision invests heavily in the field to expand its product portfolio. The present product portfolio includes state-of-the-art image enhancement software for 2D/3D/4D Ultrasound, MRI, X-Ray and Mammography which is used by leading equipment manufacturers worldwide. ContextVision is based in Sweden and listed on the Oslo Stock Exchange under the ticker COV.



**FOR MORE INFORMATION, PLEASE VISIT [WWW.CONTEXTVISION.COM](http://WWW.CONTEXTVISION.COM)**